

# Tutorial 4 – BWA, SAMtools & BCFtools Part 2

---

MICB405 – BIOINFORMATICS – 2021W-T1

08 OCTOBER 2021

AXEL HAUDUC

# BWA MEM vs BWA ALN

---

## **bwa mem syntax**

```
bwa mem ref.fasta read1.fastq read2.fastq > alignment.sam
```

## **bwa aln (“base BWA”) syntax**

```
bwa aln ref.fasta read1.fastq > read1.sai
```

```
bwa aln ref.fasta read2.fastq > read2.sai
```

```
bwa sampe ref.fasta read1.sai read2.sai read1.fastq read2.fastq >  
alignment.sam
```

# Recap of BWA ALN

---

Older, original algorithm that created SAM files

Somewhat better for read lengths below 70bp

- Tailored to older sequencers that produced  $\approx 36$ bp reads

SAI file = “suffix array index”

- Intermediate file containing Burrows-Wheeler transformed intervals of reference where each read matches
- Converted to SAM taking into account paired-end reads with secondary step “bwa sampe”

# Getting more detailed statistics on BAMs

---

`samtools stats file.bam | grep "^SN"`

samtools stats contains sets of information in lines beginning with a set of codes to denote different sections

<b>CHK</b>	Checksum
<b>SN</b>	Summary numbers
<b>FFQ</b>	First fragment qualities
<b>LFQ</b>	Last fragment qualities
<b>GCF</b>	GC content of first fragments
<b>GCL</b>	GC content of last fragments
<b>GCC</b>	ACGT content per cycle
<b>GCT</b>	ACGT content per cycle, read oriented
<b>FBC</b>	ACGT content per cycle for first fragments only
<b>FTC</b>	ACGT raw counters for first fragments
<b>LBC</b>	ACGT content per cycle for last fragments only
<b>LTC</b>	ACGT raw counters for last fragments
<b>BCC</b>	ACGT content per cycle for BC barcode
<b>CRC</b>	ACGT content per cycle for CR barcode
<b>OXC</b>	ACGT content per cycle for OX barcode
<b>RXC</b>	ACGT content per cycle for RX barcode
<b>QTQ</b>	Quality distribution for BC barcode
<b>CYQ</b>	Quality distribution for CR barcode
<b>BZQ</b>	Quality distribution for OX barcode
<b>QXQ</b>	Quality distribution for RX barcode
<b>IS</b>	Insert sizes
<b>RL</b>	Read lengths
<b>FRL</b>	Read lengths for first fragments only
<b>LRL</b>	Read lengths for last fragments only
<b>ID</b>	Indel size distribution
<b>IC</b>	Indels per cycle
<b>COV</b>	Coverage (depth) distribution
<b>GCD</b>	GC-depth

# Very quick comparison of BAMs

---

```
diff --side-by-side \  
<(samtools stats aln/bordetella.final.bam | grep "^SN") \  
<(samtools stats mem/bordetella.final.bam | grep "^SN")
```

Process substitution:

<(COMMAND) = insert output of COMMAND as a file here

- Useful when you need to pipe in 2+ items into your command

# Samtools tview

---

This is a commandline option to viewing read alignments

- Less flexible than IGV, but easy for quick comparisons

## **samtools tview syntax**

```
samtools tview file.bam --reference reference.fasta
```

# Jumping to a specific locus

---

```
samtools tview file.bam --reference reference.fasta -p  
chrom:position
```

```
samtools view file.bam chrom:start-end
```

# Filtering your BAM file

---

Select reads with a mapping quality above 40

- `samtools view file.bam -b -h -q 40 > file.filtered.bam`

Keep reads w/ matching flags (4 = UNMAPPED)

- `samtools view file.bam -b -h -f 4 > file.filtered.bam`

Discard reads w/ matching flags

- `samtools view file.bam -b -h -F 4 > file.filtered.bam`



# View alignments and vcfs together

scp reference,  
bam, bam.bai, vcf,  
and vcf.bai into  
your computer  
and load into IGV

