

Tutorial 5 – MACS2

09 OCTOBER 2020

MICB405

What is MACS2

“Peak-calling” software

Calls areas of the genome enriched for your epigenomic mark of interest based on presence of FASTQ reads

Outputs BED file which describe the genomic region where your mark is present

- This is distinct from sequence alignment!

Tutorial background

In this tutorial we will be analyzing a dataset generated as part of a project designed by Dr. Perona-Wright and Nicolette Fonesca.

These data were generated from Naive T-cells purified from mice and subjected to ChIP-seq


Specifically we will use ChIP-seq data for H3K27ac

- An active mark found at gene promoters and enhancer elements
- We will also be using an input control that represents a sequencing library generated from the sheared DNA fragments prior to immunoprecipitation

Variant Calling




Variant Calling




Let's find them!

Peak Calling



Peak Calling



These are peaks – their locations are linked to the epigenomic mark that is being studied

Conda environment

Package management environment

Allows you to install programs without having to worry about installing dependencies

We will use Conda to run MACS2

When running the downloaded shell script, remember to `chmod` to add execute privileges